

The C13366 series (GA type) are optical measurement modules capable of detecting low level light. These modules consist of a thermoelectrically cooled MPPC, a signal amplifier circuit, a high-voltage power supply circuit, and a temperature control circuit. The photosensitive area is available in two sizes of  $1.3 \times 1.3$  mm and  $3 \times 3$  mm, and the signal output is analog. Modules operate just by connecting them to an external power supply (±5 V).

## Features

- Built-in TE-cooled MPPC [MPPC for precision measurement (new product)]
- High sensitivity in the short wavelength range
- Low noise equivalent power
- Built-in temperature control function
- Analog output

## Applications

- Low-level-light measurement
- Flow cytometry
- Fluorescence measurement
- Analytical instrument

#### Structure

Parameter	Symbol	C13366-1350GA	C13366-3050GA	Unit	
Internal MPPC	-	S13362-1350DG	S13362-3050DG	-	
Effective photosensitive area	-	1.3 × 1.3	3 × 3	mm	
Pixel pitch	-	50			
Number of pixels	-	667	3600	-	

## Absolute maximum ratings

Parameter	Symbol	Condition	Value	Unit
Supply voltage	Vs		±6	V
Operating temperature	Topr	No dew condensation*1	-10 to +40	°C
Storage temperature	Tstg	No dew condensation*1	-20 to +70	°C

\*1: When there is a temperature difference between a product and the surrounding area in high humidity environment, dew condensation may occur on the product surface. Dew condensation on the product may cause deterioration in characteristics and reliability.

Note: Exceeding the absolute maximum ratings even momentarily may cause a drop in product quality. Always be sure to use the product within the absolute maximum ratings.

#### Electrical and optical characteristics (Typ. Ta=25 °C, $\lambda = \lambda p$ , Vs=±5 V, unless otherwise noted)

Parameter		Symbol	Condition	C13366-1350GA			C13366-3050GA			Unit
		Symbol		Min.	Тур.	Max.	Min.	Тур.	Max.	Unit
Spectral response range		λ		320 to 900		320 to 900			nm	
Peak sensitivity w	avelength	λр		-	500	-	-	500	-	nm
Element temperat (setting temperat		Td		-	-20	-	-	-20	-	°C
Photoelectric sense	sitivity	-		$0.7 \times 10^{9}$	$1.0 \times 10^{9}$	$1.3 \times 10^{9}$	$0.7 \times 10^{9}$	$1.0 \times 10^{9}$	$1.3 \times 10^{9}$	V/W
Cutoff froquoncy	Cutoff froquency High band	fc	-3 dB, sine wave	3	4	-	3	4	-	MHz
Cutoff frequency Low band			-5 ub, sine wave	DC			DC			-
Rise time		tr	10 % to 90%, 1.p.e.	-	5	-	-	9	-	ns
Noise equivalent	power	NEP	Dark state	-	0.1	0.2	-	0.15	0.3	fW/Hz <sup>1/2</sup>
Minimum detection	on limit	-	Dark state	-	0.25	0.5	-	0.35	0.7	pW rms
Maximum output	voltage	-		-	4.7	-	-	4.7	-	V


MPPC modules	GA type	C13366 series	
--------------	---------	---------------	--

#### Electrical characteristics

Parameter	Symbol	Condition	Min.	Тур.	Max.	Unit	
Supply voltage*2	+Vs		+4.75	+5	+5.25	V	
	-Vs		-4.75	-5	-5.25		
Current consumption	Ic	+Vs	-	+200	+1000	m۸	
		-Vs	-	-20	-40	mA	


\*2: A power supply with 1 A or higher output must be used.

# Block diagram




# Connection example

Using the supplied power cable, connect the MPPC module to a power supply. You can monitor the output waveform by connecting the MPPC module to an oscilloscope.







## Measurement example



# Photoelectric sensitivity vs. wavelength




# Frequency response (typical example)






# Linearity



# Dimensional outline (unit: mm)





#### Accessories

- · Power cable
- $\cdot$  Instruction manual

# Options (sold separately)

Coaxial converter adapter A10613 series

The A10613 series is a coaxial adapter that converts the SMB coaxial connector for signal-output on the MPPC module to a BNC or SMA coaxial connector. This adapter allows connecting a BNC or SMA cable to the MPPC module.



A10613-01 (SMB-BNC)

A10613-02 (SMB-SMA)

# Precautions

- For cleaning the product, wipe using a clean, soft, dry cloth. Do not use organic solvents such as thinner and acetone.
- Do not cover the product with a dark cloth or something similar while the product is running. Covering it can cause the internal temperature to rise and cause abnormal operation.

Type no.	Output	Effective photosensitive area (mm)	Pixel pitch (µm)	Cooling	
C13365-1350SA	Analog	1.3 × 1.3		Non-cooled	
C13365-3050SA	Analog	3 × 3			
C13366-1350GA	Angle g	1.3 × 1.3	50	TE-cooled	
C13366-3050GA	Analog	3 × 3	50	I E-COOIeu	
C13366-1350GD	Digital	1.3 × 1.3		TE-cooled	
C13366-3050GD	Digital	3 × 3			

# - Lineup of MPPC modules



# Related information

www.hamamatsu.com/sp/ssd/doc\_en.html

Precautions

Disclaimer

MPPC is a registered trademark of Hamamatsu Photonics K.K.

Information described in this material is current as of March 2020.

Product specifications are subject to change without prior notice due to improvements or other reasons. This document has been carefully prepared and the information contained is believed to be accurate. In rare cases, however, there may be inaccuracies such as text errors. Before using these products, always contact us for the delivery specification sheet to check the latest specifications.

The product warranty is valid for one year after delivery and is limited to product repair or replacement for defects discovered and reported to us within that one year period. However, even if within the warranty period we accept absolutely no liability for any loss caused by natural disasters or improper product use. Copying or reprinting the contents described in this material in whole or in part is prohibited without our prior permission.



# www.hamamatsu.com

#### HAMAMATSU PHOTONICS K.K., Solid State Division

1126-1 Ichino-cho, Higashi-ku, Hamamatsu City, 435-8558 Japan, Telephone: (81)53-434-3311, Fax: (81)53-434-5184

1120-1 ICHINO-CRO, HIGJASIN-KU, HalmarinaltSU CUty, 452-6535 Japari, Telephone: (31)53-434-3311, FAX: (31)53-434-3149 U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, N.J. 08807, U.S.A., Telephone: (1)908-231-1218, E-mail: usa@hamamatsu.com Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8, E-mail: info@hamamatsu.de France: Hamamatsu Photonics France S.A.R.L: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (49)8152-375-0, Fax: (43)169 53 71 10, E-mail: info@hamamatsu.de Northe Kingdom: Hamamatsu Photonics France S.A.R.L: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)16 95 37 10, D, Fax: (33)16 95 37 110, E-mail: info@hamamatsu.fr United Kingdom: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01, E-mail: info@hamamatsu.se Italy: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01, E-mail: info@hamamatsu.se Italy: Hamamatsu Photonics (China) Co., Ltd.: B1201, Jiaming Center, No.27 Dongsanhuan Bellu, Chaoyang District, 100020 Beijing, P.R.China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866, E-mail: hpc@hamamatsu.co.m.tm Taiwan: Hamamatsu Photonics Taiwan Co., Ltd.: 8F-3, No. 158, Section2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (86)3-659-0080, Fax: (88)3-659-0081, E-mail: info@hamamatsu.co.m.tm