Optogenetics: A Bright Future for Voltage Gated Ion Channels

AXXAA

Viviana Agus Axxam S.p.A., Milano, Italy • Hamamatsu and AXXAM: FDSS µCELL demo period

Cav1.3 and ChR2 assay

- **Optogenetics**: overview and advantages
- **Channelrhodopsin2** to modulate cell membrane voltage
- Activation of Cav1.3 by ChR2: recombinant assay setup
- Validation of "**light protocol**" at FDSS µCELL: test of reference compounds
- Comparison with "K+ protocol" and patch-clamp data
- **Conclusions** and future perspectives



FDSS µCELL DEMO @ AXXAM

Ca²⁺ assay (fluorescent dyes, luminescent photoprotein)

- CCKAR (GPCR)
- GLP1R (GPCR)
- ADORA1 (GPCR)
- DRD1-DRD2 (GPCR)
- Enzymatic assay

Glow Luminescence assay:

- PPAR α , PPAR δ (NHR)
- Promoter assay

Genetically encoded sensor

TMEM16A (EYFP)

Optogenetics

• Cav1.3 (ChR2)







Optogenetics: overview

Optogenetics is a technology that combines:

- 1) A «genetic» component, able to target specific neuron types
- 2) An **«optical»** component, able to interact specifically with the genetic component

to achieve fast control of well-defined events in specific cells of living tissue



The starting point was the idea to have a system available to control the activity of specific neuron types in the brain in a better way



Method of the year 2010



Rhodopsins: Light-gated ion channels

Opsins:

- Seven-transmembrane, light-responsive proteins
- Rapidly translocate specific ions across the membranes of the cells in which they are expressed
- Contain the Vitamin-A derived chromophore **all-trans-retinal** as a light capture molecule
- Studied since the 1970s for their fascinating biophysical properties
- Used by several different life forms that use light as energy source or sensory cue



From Zhang F. et al; Cells 147, 2011. 1446-1457

Structural simplicity, fast kinetics ► attractive tool for a rapid control of specific cellular processes, such as, for example, modulation of membrane voltage and neuronal action potentials propagation

Light as activating stimulus **b** more physiological, compared to other hyperpolarizing or depolarizing stimulus (for example K⁺ injection)

Possibility to target their expression to specific cell types **>** genetically defined modulation of cellular processes



Channelrhodopsin-2





From Ed Boyden Lab.

- Seven transmembrane opsin (eyespot of unicellular alga Chlamydomonas reinhardtii)
- Activated by blue light (470 nm): the light causes a conformational change in the light sensitive molecule (retinal), which in turn causes a conformational change and the opening of the channelrhodopsin protein
- Non-selective cation channel (Na⁺, K⁺, Ca²⁺, H⁺): the flow of ions changes the electrical potential across the cell membrane which might, if sufficiently large, cause the neuron to fire
- Widely used to depolarize neurons and generate action potential firing: very good expression in different cell hosts



From: Wong J, J Mech Phys Solids 2012 Jun 1; 60(6) 1158-1178

Optogenetic control of Cav1.3



- L-type calcium channel
- High Voltage Activated (HVA)
- $\alpha 1$ (pore) + $\alpha 2\delta$, β , γ (accessory) subunits
- Therapeutic target: Cardiovascular, hormone secretion, CNS (Parkinson's, Alzheimer's disease)
- **Drug need:** Cav1.2 selectivity; state-dependent





Cav1.3 optogenetic assay





Channelrhodopsin-2 and cell based assays

KEY QUESTIONS

- Is it possible to adapt the assay to the FDSS µCELL optics for use in HTS?
- Can ChR2 be used to depolarize cells, such as HEK293, avoiding the artificial depolarization protocols such as KCl injection?
- Does the exposure of the cells to blue light of adequate intensity induce a ChR2 dependent cellular depolarization with subsequent activation of the transfected target?
- Does the ion flux through ChR2 alter the detection of the transfected target?

GOAL

Generate stable cell lines co-expressing a Voltage Gated ion channel of interest and ChR2 without altering the ion channel pharmacology

POTENTIAL ISSUES

- The light produced by the instrument LED system might not have the adequate intensity for ChR2 activation
- The ion flux through the ChR2 might be not sufficient to induce membrane depolarization
- The membrane depolarization induced might be not sufficient to drive the activation of transfected voltage gated channels.
- ChR2 is not permeable to Ca²⁺ in the presence of extracellular Na⁺; therefore Cav channels are ideal targets to be modulated with optogenetics, since their activity can be monitored by the use of a Ca²⁺ sensitive dye



ChR2 induced membrane depolarization





Cav1.3 half-inactivation protocol





State-dependent blockers with «Light protocol»



Light protocol vs. K⁺ protocol vs. qPatch





Summary and conclusions

Isradipine IC ₅₀	"K⁺ protocol"	"Light protocol"	qPatch 16x	Literature
Resting	95 nM	137 nM	362 nM	300 nM (-90m∨)
Half-inactivated	14 nM	12 nM	32 nM	30 nM (-50m∨)
RATIO	6.8	11.4	11.3	10

MAIN ACHIEVEMENTS:

- FDSS µCELL optics is well suitable for ChR2 activation
- A "Light protocol" was set up at the FDSS µCELL to study the Cav1.3 channel either in resting or inactivated state
- The pharmacology of known state dependent blockers has been successfully validated, showing a good agreement with the classical "K⁺ protocol", patch clamp experiments and literature data

HIGHLIGHTS:

- First time ChR2 used for optical control of recombinant voltage-gated calcium channel assay
- Physiological, robust, precise activation of Cav1.3 channel

FUTURE PERSPECTIVES:

• Light modulation of other voltage-gated ion channel target is ongoing





Aknowledgments

AXXAM:

Alberto di Silvio cell line generation Sara Tremolada cell line validation Jean-Francois Rolland patch-clamp Katharina Montag clonings Loredana Redaelli cell biology head Lia Scarabottolo discovery services director

Stefan Lohmer overall strategies

<u>Hamamatsu team:</u> Jean Marc d'Angelo Annamaria Mauro Laura Confalonieri

Via Meucci 3 20091, Bresso (Milan, Italy) phone + 39 02 210561 fax + 39 02 2105602 www.axxam.com