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Introduction

Please read carefully

- Please carefully read the “Precautions” before using.
- Please read the attached manual carefully and use the product correctly.
- Please keep the attached manual and read when it is necessary.

This document uses the following safety alert symbols and signal words.

 DANGER	<p>Indicates an imminently hazardous situation which if not avoided, will result in death or serious injury. This signal word is to be limited to the most extreme situations.</p>
 WARNING	<p>Indicates a potentially hazardous situation which if not avoided, could result in death or serious injury.</p>
 CAUTION	<p>Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It is permitted for property damage only accidents.</p>

 Note	<p>Indicates a note to help you get the best performance from the system. Read the contents of the note carefully to ensure correct and safe use. Failure to observe one of these notes might impair the performance of the system.</p>
	<p>Indicates a cautionary item that should be obeyed when handling the system. Read the contents carefully to ensure correct and safe use.</p>
	<p>Indicates an action that is forbidden. Read the contents carefully and be sure to obey them.</p>
	<p>Indicates a compulsory action or instruction. Read the contents carefully and be sure to obey them.</p>

Safety consideration

Make sure you follow the instructions

This product is intended for use in general electronic equipment (Measurement equipment, office equipment, information and communications equipment, home appliances, etc.). Do not use this product for special purposes (Nuclear control equipment, Aerospace equipment, Medical equipment and Transport equipment that directly affect human life, disaster prevention and safety equipment, etc.) that require extremely high reliability and safety, except as described in the individual product documentation. In order to use this product correctly and safely, make sure to follow the safety notices below when operating this product. Hamamatsu Photonics is neither responsible nor liable for damages caused by using this product contrary to the safety notices.

WARNING

● In case an abnormality has occurred

Do not use this product if smoke, unusual odor, abnormal noise, or heat is generated during operation. Continuing operation may cause a fire or electrical shock. In this case, immediately turn off the PC (personal computer) connected to this product, and then unplug the USB cable. After making sure that the abnormal condition has disappeared, contact our sales office.

Never attempt to repair this product yourself since it can be dangerous.

● Do not use this product above its absolute maximum ratings

Do not use this product above values stated in the absolute maximum ratings. Using this product beyond its absolute maximum ratings may cause damage or fires.

● Absolutely do not make mistaken connections

Doing so may cause overcurrent, causing damage or fires in parts and devices.

● Do not damage cords and cables

Do not scratch, damage, modify, or place heavy objects on top of cords. Doing so can break cords and cause electric shock or fire.

● Make sure to handle the power cable

Make sure to correctly connect the power cable to this product and the power cable to an external power supply. Improper handling or connection may cause electric shock or fire, which is extremely dangerous.

● Make sure to connect a compatible cable

Make sure to use a compatible cable when using this product. This product is designed to connect with a power cable using LEMO connector FGG.00.302.CLAD35. Using an inappropriate connector or inappropriate cable may cause electric shock or fire, which is extremely dangerous.

- **Make sure to install in the correct location**

Make sure to install this product in the correct location, such as away from dust, away from splash water, and in a well-ventilated area, to prevent fire or power failure.

- **Do not use in gas**

Do not use this product in the presence of flammable or explosive gases or vapors. Using this product in such an environment is extremely dangerous.

- **Do not vibrate or impact**

Excessive vibrations or impacts may cause damage or misalignment to parts, resulting in fire, electric shock, or failure.

- **Do not enter foreign objects or water**

It may cause fire or electric shock if flammable objects, metals, or moisture enter this product, which is extremely dangerous.

- **Do not disassemble**

Do not disassemble this product. Adjusting or modifying the equipment may cause fire, electric shock, burns, or malfunction, etc.

- **Do not look directly at the emitted light**

This product uses a tungsten lamp system. Light in the wavelength range visible to near-infrared is emitted. Emitted light is classified as exempt in the light risk group: EN 62471 / IEC 62471, but it does not guarantee that there will be no impairment of visual function when viewed directly or at close range. Do not look directly at the emitted light.

- **Beware of high temperature**

The temperature of this product rises after the lamp on. If you touch this product directly for a long time after the lamp on, it may cause redness, itching, rashes, etc. on your skin, or it may cause burns. Before handling this product, turn off the lamp and make sure this product is sufficiently cooled.

■ Notices on handling this product

⚠ CAUTION

In order to use this product safely, Users should prepare the necessary peripheral equipment and follow the precautions regarding the handling of the products.

Use this product within the absolute maximum ratings.


Exceeding the absolute maximum rating even momentarily may degrade the quality of this product. Make sure to use it within the absolute maximum ratings. When applying voltage to this product, be careful not to exceed the absolute maximum rating.


Do not apply overvoltage to this product.

If a voltage exceeding the absolute maximum rating (6V) is applied to this product even momentarily, the filament of the built-in lamps will melt, deform, or short circuit, then it will cause the lamp not to light. Also set the limit of the absolute maximum rating (550mA) on the applied current to prevent overvoltage. If damage due to overvoltage is confirmed, repair will be charged even if it is within the product warranty period. Please contact us for lamps replacement.

Example of normal lamp filament

Example of damaged filament due to apply overvoltage (lamp does not light)

Use a suitable external power supply.

Use a DC stabilized power supply that the light intensity can be adjusted according to the measurement samples or measurement environments. Make sure to connect this product, the power cable, and the external power supply correctly.

Attach a heatsink to this product to exhaust heat.

After the lamp on, the temperature of this product gradually rise, reaching around 60°C, when this product is used alone at room temperature. In order to suppress the temperature rise, make sure to use this product with a heatsink or cooling fan. It is recommended to insert a material with high thermal conductivity (grease, sheet, etc.) between this product and the heatsink.


Tighten this product correctly with screws.

There are two mounting holes for M3 screws on the rear panel of this product. Make sure to tighten with the specified torque for M3 screws. The specified screw torque is set according to the material of the receiving side. Also, tighten the two screws alternately that a uniform force is applied. Inappropriate screwing may cause problems such as fixing this product.


Connect the cable with moderate force.

Applying excessive force or any force in other direction may cause poor connection or damage to the connector. The connector of this product is EGG.00.302.CLL (LEMO). It is a two-core connector. Also make sure to use the connector: FGG.00.302.CLAD35 (LEMO) for the cable. For connection, use the red line as a mark, and inset all the way in. Make sure to refer to the LEMO catalog for details.

Plug in

Plug in

completed

Use the optical fiber with handling precautions.

Clean the tip of the optical fiber before connecting this product. Absolutely avoid inserting a protruding object into the optical fiber connector. It may damage the optical fibers. Also, use a compatible optical fiber for connecting this product. If you use a non-compliant optical fiber, you may not be able to obtain the full characteristics of this product, or it may cause problems of this product.

Store this product in a clean environment.

Do not store or operate this product in the environment such as where is a dusty, watery, or condensation. Also, store this product with the protective cap which is attached at the time of shipment.

Do not touch the light exit port with bare hands.

The light output intensity decreases due to devitrification such as burning of dirt, if lamp is on with dust or fingerprints on it. If the lamp is dirty with finger marks, etc., wipe it off with gauze or absorbent cotton that has been soaked in ethyl alcohol and tightly wrung out, before turning on the lamp.

Do not apply excessive vibration, impact, or force.

Excessive vibration, shock, or force can cause such as misalignment of internal optics, damage to the window, damage to the internal lamps, break to the filament of lamps, or failure to the lamps.

Not designed for use in high places.


This device is not designed for use in high elevation places. Do not use it in places where the elevation is over 2,000 meters.

Please contact our help desk to replace the lamp.

Using the lamp for a long time will consume the filament of the lamp and eventually cause it to burn out. Hamamatsu Photonics will do alignment and electrical connection required for lamp replacement, so please contact us to do replacement or repair of the lamp in this device. In order to obtain the full characteristics of this product, if even one of the four lamps in this product does not light up, please contact us for a repair request.

Four lamps on


One lamp does not light

The lighting can be confirmed in the emission state.

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