


Digital CMOS board level camera C13949-50U



The core of the C13949-50U camera is the scientific image sensor, an advanced CMOS detector that simultaneously achieves high resolution, fast readout speeds and low noise. With its small size, board-level design and simple, low-cost integration using USB 3.0, the C13949-50U is the ideal camera for OEM scientific imaging.

It is available to offer customization following your request. Please contact us for details.


FEATURES

- Resolution: 12 megapixels (4096 x 3008)
- Readout noise: 2.3 electrons (Typ.)
- Frame rate: 15 frames/s (Full resolution)
- Global shutter


APPLICATIONS

- Fluorescence imaging
- DNA sequencing
- CTC detection
- X-ray I.I. readout
- X-ray CT

SPECTRAL RESPONSE



SYSTEM CONFIGURATION






SPECIFICATIONS

Imaging device Effective number of pixel		
		CMOS image sensor
	ls	4096 (H) × 3008 (V)
Pixel size		3.45 μm (H) × 3.45 μm (V)
Effective area		14.13 mm (H) × 10.37 mm (V)
Full well capacity (Typ.)		10 500 electrons
Analog gain		1× to 15.8× (240 steps)
Dynamic range (Typ.) *1		4565 : 1
Frame rate (Internal trigger mode)	Full resolution	15 frames/s
	2×2 Binning	57 frames/s
	2×2 Binning, 4 Line Sub-array	1540 frames/s
Readout noise (rms Typ.)		2.3 electrons (Analog gain 1×)
Digital output		12 bit
Exposure time	Internal trigger mode	25.45 μs to 1 s
	External trigger mode	25.45 µs to 1 s (Edge trigger, Start trigger) 1 ms to 1 s (Level trigger, Synchronous readout trigger)
External trigger input mode		Edge trigger, Level trigger
		Synchronous readout trigger, Start trigger
Trigger delay function		0 s to 10 s in 1 µs steps
Trigger output		Programmable timing outputs
		Trigger ready output
		Global exposure timing output
		Low output, High output
Binning *2		2×2
Sub-array Sub-array		Yes
Interface		USB 3.0
Lens mount		C-mount
Connector		USB 3.0 microB type
Power supply		DC +12 V
Power consumption		Approx. 6.5 W
Ambient operating temperature		0 °C to +40 °C
Ambient operating humidity		30 % to 80 % (with no condensation)
Ambient storage temperature		-10 °C to +50 °C
Ambient storage humidity		90 % max. (with no condensation)

^{*1:} Calculated from the ratio of the full well capacity and the readout noise.

DIMENSIONAL OUTLINES

(Unit: mm)



OPTIONS

- Board camera case A12599-03
- External trigger cable SMA-BNC 5 m A12106-05
- External trigger cable SMA-SMA 5 m A12107-05
- External trigger cable 5 m A11892-05

- Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.
- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.

© 2020 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-433-8031, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH.: Arzbergerstr. 10, D-82211 Herrschling am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy, 04882 Massy, 04882

North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.set Italy: Hamamatsu Photonics (China) Co., Ldt.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Beliu, China: Hamamatsu Photonics (Taila: China: Hamamatsu Photonics (China) Co., Ldt.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Beliu, Choyang District, (10020 Beijing, PR. China: Helphone: (86)10-6586-8006, Fax: (86)10-6586-8066 E-mail: hpc@hamamatsu.com.cn
Taiwan: Hamamatsu Photonics Taiwan Co., Ltd.: 8F-3, No.158, Section 2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

^{*2:} Vertical 2 pixels are added inside the CMOS image sensor, and horizontal 2 pixels are added digitally inside the camera.